HANDLEIDING SUBSTRATES & SUBSTRATE BUFFER		
product code: R007R.0000 R008R.0000 R010R.0000		
SNUJDEN EN PLAKKEN VAN PARAFFINE COUPES PROTEOLYTISCHE VOORBEHANDELING HYBRIDISATIE PROCEDURE		

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**Technical specifications**

Cat. No.	Description	Contents	Format	Appearance	pH	Result
R007R.0000	3-amino-9-ethylcarbazole HRP substrate for use in conjunction with AEC buffer (SUBS AEC)*	2 mL (100 tests), clear vial, blue cap	Concentrated in dimethylformamide	Dark brown liquid	NA	Reddish-brown precipitate
R010R.0000	AEC buffer solution containing hydrogen peroxide for use in conjunction with AEC substrate (BUFF AEC)*	15 mL (100 tests), white vial, blue cap	Pre-diluted	Clear liquid: colourless	4.8 ± 0.2	Reddish-brown precipitate
R008R.0000	Nitroblue tetrazolium salt and 5-bromo-4-chloro-3-indolyl phosphate AP substrate (SUBS NBT/BCIP)	15 mL (100 tests), white vial, blue cap	Ready to use	Clear liquid, yellow	9.0 ± 0.2	Indigo precipitate

\*Directions for preparing AEC working solution:

# of specimens	# of drops of AEC substrate	Vol. of AEC buffer
1-13	4	2 mL
14-26	8	4 mL
27-39	12	6 mL
40-52	16	8 mL

Application : substrate for colorimetric detection of conjugates by *in situ* hybridisation (ISH)  
 Performance : 3.2 pg DNA (labelled and conjugated), checked by direct spot blot  
 Storage : refrigerated (2-8 °C) protected from light  
 Stability : until expiry date printed on label  
 Precautions : - AEC is a potential carcinogen  
                   - homogenise solutions before use  
                   - avoid contact with eyes and skin; do not swallow

**Related products**

Please contact your local supplier for further information.

Purchase does not include the right to exploit this product commercially and any commercial use without the explicit authorization of PanPath BV is prohibited.

## Limitations of Procedure

### Product Precipitating Substrates

- The REMBRANDT® Precipitating Substrates are solely applicable for the visualisation of conjugated antibodies which are used for the detection of hybridised Biotin or Digoxigenin labelled DNA and RNA probes in ISH, which may be present in cell preparations (paraffin sections, frozen sections or cytological specimen).
- Appropriate medical decisions are only possible if the medical traceability is ensured. The product is intended for professional use as an aid in the diagnosis by *in situ* hybridisation.
- Either human tissue sections or human cytological preparations may be used. Samples must be fixed in buffered formalin or alcohol. Sections should be cut at 4 µm thickness, glued to the glass slides with a bio-adhesive (e.g. organosilane), dried at room temperature, subsequently dried at 37 °C overnight, complete deparaffinisation in xylene and alcohol series and air dried. Cytological specimen should be prepared as required by the user, fixed with cytological fixation agent, rinsed in distilled water prior to the ISH procedure and air dried.
- Many factors can influence the performance of the ISH procedure. Failure in detection can be due to i.e. improper sampling, handling, the time lapse between tissue removal and fixation, the size of the tissue specimen in the fixation medium, the fixation time, processing fixed tissue, the thickness of the section, the bio-adhesive on the slide, deparaffinisation procedure, incubation times, incubation temperatures, all other reagents (i.e. deparaffinisation reagents, proteolytic treatment, probes etc.) used in the procedure and interpretation of results.
- The performance of the ISH procedure is also affected by the sensitivity of the method and the DNA or RNA target load; in case the limit of the sensitivity is reached or when the target DNA or RNA load is too low, a false negative reaction may be the result.
- The clinical interpretation of the results should not be established on the basis of a single test result. Moreover, diagnosis should also take the clinical history, symptoms, as well as morphological data into consideration. Negative results therefore do not rule out any possibility of a positive specimen.
- The REMBRANDT® test results are not to be relied on in case the sampling, sampling method, quality, sample preparation, other reagents used, controls and procedure followed is not optimal.
- Therapeutic considerations based on the result of this test alone should not be taken. Positive results should be verified by other traditional diagnostic methods such as but not limited to clinical history, symptoms, as well as morphological data.
- The medical profession should be aware of risks and factors influencing the intensity, the absence or presence of probe signals which can not be foreseen when applying this test.
- The user should carefully consider the risk and use of sample material for this test in case the sample material does not contain sufficient or representative test material.
- Laboratory personnel performing the test should be knowledgeable and be able to interpret the test results.

### Interpretation of the results

First, check the negative and positive controls that have been incubated with the test slides simultaneously:

- The negative control should be really negative, i.e. not show any localised colour precipitations. If the negative control could be interpreted as being positive, discard the results since no conclusions can be drawn.
- The positive control should show colour precipitations in conformity with the localisation of the target DNA or RNA. The colour should show the proper shade and must be clearly visible in the preferential cell/ tissue type and correspond to the target localisation.

In the test slides, start under low power magnification and focus on localisation and colour to see whether:

- The positivity (colour precipitation) observed is localised in the cell type preferred by the target.
- The colour has the right shade (no endogenous or formalin pigment).

Use high power magnification to see whether:

- The positive staining texture (granular, etc), demarcation and localisation are conform the positive control staining pattern.

### Product in combination with other devices

The REMBRANDT® Precipitating Substrates are intended for stand-alone usage. The in vitro diagnostic is intended to be used in combination with standard formalin fixed, paraffin embedded tissue blocks, standard tissue freezing, tissue sectioning (microtome), standard cytological preparation methods, hot plate(s), stove(s), incubation device(s), water bath(s), temperature and incubation time control(s), other needed reagents (but not supplied with this reagent) for ISH and a microscope. The combination has been tested and validated. Since the standard formalin fixed, paraffin embedded tissue blocks, standard tissue freezing, tissue sectioning (microtome), standard cytological preparation methods, hot plate(s), stove(s), incubation device(s), water bath(s), temperature controls, incubation time control(s) and other needed reagents such as but not limited to proteolytic reagents, labelled probes, conjugates and a microscope is not combined with the device as a product, conformity with the essential requirements is not applicable. Assay validation criteria are mentioned in '*Interpretation of the results*' and are also depending on the target; since the Precipitating Substrates can be used for all targets labelled with AP or HRP conjugates and the results also depend on the target load , the validation criteria may vary.

### References

- Graham, R.C., et al., J. Histochem. Cytochem., 13, 150 (1965).
- Kaplow, L.S., am J. Clin. Path., 63, 451 (1974).
- Blake, M., Anal. Biochem., 136, 175 (1984).
- Horowitz, J., et al., J. Med. Chem., 9, 447 (1966).

## REFERENCE GUIDE SUBSTRATES & SUBSTRATE BUFFER

product code: R007R.0000  
R008R.0000  
R010R.0000

### PRETREATMENT OF PARAFFIN SECTIONS PROTEOLYTIC TREATMENT HYBRIDIZATION PROCEDURE

Start with *DETECTION PROCEDURE* (conjugate) and proceed with A or B:

#### (A) AEC substrate

- Prepare AEC (blue) work solution according the following table
- | # of specimens | # of drops of AEC substrate | Vol. of AEC buffer |
|----------------|-----------------------------|--------------------|
| 1-13           | 4                           | 2 mL               |
| 14-26          | 8                           | 4 mL               |
| 27-39          | 12                          | 6 mL               |
| 40-52          | 16                          | 8 mL               |

- Apply 2-3 drops of the work solution to each specimen and incubate in dark
- Tap off excess substrate solution and wash slides in distilled/deionised water

#### (B) NBT/BCIP substrate

- Apply 2-3 drops of NBT/BCIP substrate(blue) to each specimen and incubate in dark
- Apply 2-3 drops of the work solution to each specimen and incubate in dark
- Tap off excess substrate solution and wash slides in distilled/deionised water

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